



Weight of frame 30 kg - 66 lbs.
Weight of unit 80 kg - 176 lbs.
Total weight 110 kg - 242 lbs.



Note: See DPA-2522-10

Note: See DPA-2522-10

Sanitary Installation	MEIKO	In the building
B Cold Water	R 1/2" NPT, Stop-cocks, 1/4" Strainer	Inlet pipe 3/4", 41"-77"
C Hot water	R 1/2" NPT, Stop-cocks, 1/4" Strainer	Inlet pipe 3/4", Stop-cocks 110"-140"
E Waste Outlet	4" TRAP	Inlet pipe 3/4", Stop-cocks 24.76 gal./minute, 14.45 PSI
Electrical Installation		
A 15 Ph - 208/230V, 60Hz, 4.2kW/1.6kVA		
Protection class Equipment class 1, Pollution degree 2		
Protection class according to IEC 384-5-34 and conform to the local and national power utilities regulations		
Connection dimensions are subject to tolerances. It is neither allowed to give this drawing without our permission nor may it be used for competition purposes. All rights are reserved. Subject to technical changes.		
Reference	MEIKO USA INC. 1349 Hill Quaker Blvd. LaVerne TN 37086 TEL: (615)399-6600 FAX: (615)399-6600	
MD Anderson Cancer Center	Altert Drawn	Checked
Type TDP LINE 30 DT -NDV	Thermal disinfection/intermittent steam generation	
Proj. no.		
Drawing no.		
Scale: 1:60	Drawn: cm	Date: Dc-30.2007
	Checked: cm	Date: Dc-30.2007
		Index